



For more information, please contact:

info@chrom4.com | www.chrom4.com

Cyclic Hydrocarbons

Column: **CH1 - 10m, 0.53mm, 5.00 μ m**

Catalog Code: **CH1-1053500**

Conditions:

Injection: Split 250°C, 1.0 μ L.

Detector: FID 250°C.

Oven Temp.: 40°C, 5°C/min, 100°C.

Carrier Gas: Helium, 5mL/min.

Peaks:

1. Cyclohexane
2. Cycloheptane
3. Cyclooctane
4. n-Decane

