



For more information, please contact:

info@chrom4.com | www.chrom4.com

Drugs of Abuse - Basic Drugs

Column: **CH1 - 15m, 0.25mm, 0.25 μ m,**

Catalog Code: CH1-152525

Conditions:

Injection: Split 280°C, 1.0 μ L, 1:50 Split Ratio.

Detector: FID 320°C.

Oven Temp.: 120°C, 25°C/min, 310°C.

Carrier Gas: Hydrogen, 30 kPa.

Peaks:

1. Caffeine
2. Lidocaine
3. Procaine
4. Cocaine
5. Butyl Antrachinon (I.S.)
6. Heroine
7. Papaverine
8. Etaverine (I.S.)
9. Narcotine